


Ringbio




Companion Animal Test Kit

Rapid Test Kit, RT-PCR ELISA, FIA Tests

rapid and easy-to-use kit to monitor your beloved pets
for veterinary clinic and home use **

Ring Biotechnology Co Ltd


** may not apply in certain countries



Lateral flow rapid test kits


Item #	Product	Description
C11281	Canine parvovirus / CPV antigen test card	15 Tests / box
C11291	Canine distemper virus / CDV antigen test card	15 Tests / box
C11371	Canine coronavirus / CCoV antigen test card	15 Tests / box
C11421	Canine adenovirus 2 / CAV-2 antigen test card	15 Tests / box
C11431	Canine pancreatitis test card	15 Tests / box
C11441	Canine progesterone test card	15 Tests / box
C11451	Feline Herpes Virus / FHV Antigen test card	15 Tests / box
C11361	Feline panleukopenia virus / FPV antigen test card	15 Tests / box
C11391	Feline Caliciviurs / FCV antigen test card	15 Tests / box

Test result in 10min after simple sample extration




You can always visit www.ringbio.com for our latest product update
 Not all test kits are registered in all countries, consult your sales before order
 For veterinary use only, for research use only, for outside of China only.

Fluorescence Rapid Test kits



Ringbio



Fluorescence Rapid Test kits


Item #	Product	Description
F10311	Canine C Reactive Protein (CRP) FIA Test Card	15 Tests
F10461	Feline Serum Amyloid A (fSAA) FIA Test Card	15 Tests
F11281	Canine Parvovirus (CPV) Antigen FIA Test Card	15 Tests
F11291	Canine Distemper Virus (CDV) Antigen FIA Test Card	15 Tests
F11292	Canine Coronavirus Virus (CCoV) Antigen FIA Test Card	15 Tests
F11361	Feline Panleukopenia Virus (FPV) Antigen FIA Test Card	15 Tests
F11362	Feline Coronavirus Virus (FCoV) Antigen FIA Test Card	15 Tests
F11361	Feline Herpes Virus (FHV) Antigen FIA Test Card	15 Tests
F11363	Feline (FPV/FHV/FCV) 3X antibody FIA Test Card	15 Tests
F11364	Canine (CHV/CPV/CDV) 3X antibody FIA Test Card	15 Tests
F11365	Helicobacter pylori Antigen Fluorescence Test Kit	15 Tests
F1001	FIA Strip Reader	1 set

**All kits above are based on fluorescence immunoassay method.
The FIA Reader is necessary for result determination


You can always visit www.ringbio.com for our latest product update
Not all test kits are registered in all countries, consult your sales before order
For veterinary use only, for research use only, for outside of China only.

Ring Biotechnology Co Ltd

ELISA Kits



Ringbio



ELISA Kits

Item #	Product	Description
KE0011	Canine Rabies IgG Antibody ELISA Kit	96 Tests
KE0012	Feline Rabies IgG Antibody ELISA Kit	96 Tests
KE0013	Canine Hep/Parvo/CDV Antibody Triplex ELISA Kit	96 Tests
KE0014	Feline FPV/FHV/FCV Antibody Triplex ELISA Kit	96 Tests
KE0015	Canine Hep IgG Antibody ELISA Kit	96 Tests
KE0016	Canine Parvovirus IgG Antibody ELISA Kit	96 Tests
KE0017	Canine Distemper Virus, CDV Antibody ELISA Kit	96 Tests
KE0018	Feline FPV IgG Antibody ELISA Kit	96 Tests
KE0019	Feline FHV IgG Antibody ELISA Kit	96 Tests
KE0020	Feline FCV IgG Antibody ELISA Kit	96 Tests


** All kits above are based on enzyme linked immunoassay method.
ELISA Reader is necessary for result determination

You can always visit www.ringbio.com for our latest product update
Not all test kits are registered in all countries, consult your sales before order
For veterinary use only, for research use only, for outside of China only.

Real time PCR kits



Ringbio



Real time PCR kits

Item #	Product	Description
PZ1001	Feline Panleukopenia Virus FPV Real-time PCR	16 Tests
PZ1002	Feline Coronavirus Real-time PCR	16 Tests
PZ1003	Feline Herpesvirus Real-time PCR	16 Tests
PZ1004	Feline Calicivirus Real-time PCR	16 Tests
PZ1005	Chlamydophilia felis Real-time PCR	16 Tests
PZ1006	Mycoplasma felis Real-time PCR	16 Tests
PZ1007	Feline infectious peritonitis virus Real-time PCR	16 Tests
PZ1008	Feline Chlamydophilia felis & Mycoplasma felis Real-time PCR	16 Tests
PZ1009	Feline Coronavirus & Infectious peritonitis Real-time PCR	16 Tests
PZ1010	Feline Calicivirus & Herpesvirus Real-time PCR	16 Tests
PZ1011	FCV & FHV & C. felis & M. felis Real-time PCR	16 Tests

** All PCR assay kit above are one step test kit without extraction step needed

You can always visit www.ringbio.com for our latest product update
Not all test kits are registered in all countries, consult your sales before order
For veterinary use only, for research use only, for outside of China only.

Ring Biotechnology Co Ltd

Real time PCR kits



Ringbio

Item #	Product	Description
PZ1012	Canine Distemper Virus CDV Real-time PCR	16 Tests
PZ1013	Toxoplasma gondii Real-time PCR	16 Tests
PZ1014	Babesia canis Real-time PCR	16 Tests
PZ1015	Bartonella henselae Real-time PCR	16 Tests
PZ1016	Feline Leukemia Virus Real-time PCR	16 Tests
PZ1017	Feline Immunodeficiency Virus Real-time PCR	16 Tests
PZ1018	Feline Hemoplasmas(Classify) Real-time PCR	16 Tests
PZ1019	Mycoplasma canis Real-time PCR	16 Tests
PZ1020	Canine Parvovirus Real-time PCR	16 Tests
PZ1021	Canine Adenovirus type-1 Real-time PCR	16 Tests
PZ1022	Canine Adenovirus type-2 Real-time PCR	16 Tests
PZ1023	Canine Infuenza Virus Real-time PCR	16 Tests
PZ1024	Canine Parainfluenza Virus Real-time PCR	16 Tests
PZ1025	Bordetella bronchiseptica Real-time PCR	16 Tests
PZ1026	Leptospira Real-time PCR	16 Tests
PZ1027	Tritrichomonas foetus Real-time PCR	16 Tests
PZ1028	Giardia lamblia Real-time PCR	16 Tests
PZ1029	Canine Rabies Virus Real-time PCR	16 Tests
PZ1030	Canine Coronavirus Real-time PCR	16 Tests
PZ1031	Feline Hemoplasmas Universal Real-time PCR	16 Tests
PZ1032	FeIV & Leukemia & Bartonella henselae Real-time PCR	16 Tests
PZ1038	Babesia Canis and Babesia Gibsoni Real-time PCR	16 Tests
PZ1039	Canine CoV & CPV & Giardia & Tritrichomonas foetus PCR	16 Tests
PZ1040	Canine CCoV & CPV Real-time PCR	16 Tests
PZ1041	CDV & CPIV & CAV2 & Bordetella bronchiseptica PCR	16 Tests
PZ1042	TLGRBH 6in1 Real-time PCR	16 Tests
PZ1043	FeIV & FeLV & Bartonella henselae & Hemoplasmas PCR	16 Tests
PZ1047	Canine Richettsia Real-time PCR Kit	16 Tests

** All PCR assay kit above are based on real-time fluorescence qPCR.

You can always visit www.ringbio.com for our latest product update
Not all test kits are registered in all countries, consult your sales before order
For veterinary use only, for research use only, for outside of China only.

Ring Biotechnology Co Ltd



Ringbio

NBGen Q1600 Real time fluorescence PCR Cycler



Portable
Quantitative
Accurate


Main Features

- Dual channel design with 8 wells per each
- Suitable for 2 different assays
- Functional software for quantitative analysis
- Portable with DC/AC power support
- 7inch touch screen
- 20Gb flash drive for more info storage
- WIFI/USB connection

Sample tube	0.2ml, 2 x 8 tubes,
Reaction volume	10 - 100 ul
Temp. range	0-100°C
Heating speed	6 °C /s
Temp. precision	0.1°C
Temp. control	±0. 25°C
Wave length	F1 520-540nm / F2 540-580nm
Dyes applied	F1 FAM/SYBR Green F2 HEX/VIC/JOE/TET
Light source	LED
Detection range	0 - 10 ¹⁰ copies
Dimension	300x267x198mm(L x W x H)
Weight	5.0 KG


You can always visit www.ringbio.com for our latest product update
Not all test kits are registered in all countries, consult your sales before order
For veterinary use only, for research use only, for outside of China only.

Ring Biotechnology Co Ltd



Ringbio

NBGen K3 Plus ELISA Reader



Portable
Quantitative
Accurate

7 inch colored touch screen

Superior flexibility and reliability

Excellent optical system

Extended memory measurement program

Practical and simple software

Specifications

Parameter	Details	Parameter	Details
Dimension(W*H*L)	140*420*320mm	OD precision	0.0001 A
Weight	11 Kg	Linear range	OD: 0 ~ 3.00
Power	220V, 50Hz	Linearity	R ² >0.995
Transportation	-40°C - 70 °C	Accuracy	0-0.3A, +/-0.007
Storage	-25°C - 50 °C	Accuracy	0.3-2.0A, +/-1%
Plate type	96 wells plate	Precision	C.V.<1%
Optical system	QHT lamp(8v/50w)	Measure time	8s
Weave length	400-750nm	Measure time	20s in step mode
Filter	405/450/492/630nm	Shaking	3 modes
Optional filter	8 filters in total	UI screen	7' touchscreen
Measure range	0 ~ 4.0A	Operation	by PC or touch
Filter half bandwidth	3 - 9nm	Storage	500 programs

Ringbio

-  **Ring Biotechnology Co Ltd**
-  **Building 3, Zhongtongtai TechnoPark, BDA, Beijing 100176, China**
-  **+86-10-56267496 / +86-10-26242030**
-  **info@ringbio.com**
-  **www.ringbio.com**