

Ringbio

ONE STEP RAPID TEST KIT



heart of love



*One bottle of milk
One bottle of health*



Camel milk fraud test kit

Camel milk test kit


to detect goat/sheep & cow milk in camel milk

Product code: 420003
Unit size: 96 Tests

Camel milk test kit

to detect cow milk in camel milk

Product code: 410003
Unit size: 96 Tests



Goat / sheep milk fraud test kit

Cow milk test kit

to detect cow milk content in goat/sheep milk

Product code: 400003
Unit size: 96 Tests

Ringbio is a professional manufacturer of rapid tests for food and animal diagnostics. The company is ISO9001, ISO13485 and ISO 14001 certified with GMP licensed for diagnostic kit research and production.

Ringbio

Ring Biotechnology Co Ltd



Camel milk fraud test kit

Camel milk test kit

to detect goat/sheep & cow milk in camel milk

Product code: 420003

Unit size: 96 Tests

Kit components:

Reagent Tube, 12 pcs, In each tube, 8 test strips & 8 microwell reagent

Micropipette, 100ul, 1 set

Yellow micropipette tips, 100pcs

Kit instruction manual, 1 set

Limit of detection: Cow milk 0.3% - 0.5% by volume

Goat / sheep milk 0.1 - 0.2% by volume

Camel milk test kit

to detect cow milk in camel milk

Product code: 410003

Unit size: 96 Tests

Kit components:

Reagent Tube, 12 pcs, In each tube, 8 test strips & 8 microwell reagent

Micropipette, 200ul, 1 set

Yellow micropipette tips, 100pcs

Kit instruction manual, 1 set

Limit of detection: Cow milk 0.2% - 0.3% by volume

Goat / sheep milk fraud test kit

Cow milk test kit

to detect cow milk content in goat/sheep milk

Product code: 400003

Unit size: 96 Tests

Kit components:

Reagent Tube, 12 pcs, In each tube, 8 test strips & 8 microwell reagent

Micropipette, 200ul, 1 set

Yellow micropipette tips, 100pcs

Kit instruction manual, 1 set

Limit of detection: Cow milk 0.2%-0.3% by volume



Scan the QR code to view the training video

Visit our Youtube Channel <https://www.youtube.com/c/ringbio>